

Carriers for Medium Optical Rails Feet for Large Optical Rails

CAA-LS/MB
OBT


CAA-LS/MB

RoHS Catalog Code W6024


Carriers mount on a medium optical rail, slide freely and can be locked in place at any location.

Outline Drawing


CAA-25LS/25LSEE/25LSUU




CAA-40LS




CAA-60LS




CAA-65LSEE/65LSUU



CAA-100LS



CAA-60MB



Specifications

Part Number			Compatible Optical Rails	Weight [kg]
M4	M6	INCH	Medium optical rails (OBA-SH)	0.1
CAA-25LS	CAA-25LSEE	CAA-25LSUU	Medium low-profile optical rails (OBT-SH/SUU)	0.13
CAA-40LS	-	-	Reference D024	0.17
CAA-60LS	CAA-65LSEE	CAA-65LSUU		0.19
CAA-100LS	-	-		0.28
CAA-60MB	-	-		0.18

Primary material: Aluminum
Finish: Black Anodized

OBT


RoHS Catalog Code W6026




● OBT is a thin plate type of Rail Feet; it is made to be matched with OBT-LH/SH optical rail. It allows +/-10mm height adjustment.

Outline Drawing


OBT-LS



OBT-LW



OBT-WL



Specifications

Part Number	Compatible Optical Rails	Weight [kg]
OBT-LS	Large low-profile optical rails (OBT-LH) Reference D022	0.28
OBT-LW	Medium low-profile optical rails (OBT-SH) Reference D024	0.33
OBT-WL	Large low-profile optical rails (OBT-LH) Reference D022	0.24
	Medium low-profile optical rails (OBT-SH) Reference D024	

Primary material: Aluminum
Finish: Black Anodized

Application Systems

Optics & Optical Coatings

Holders

Bases

Manual Stages

Actuators

MotORIZED Stages

Light Sources

Index

Guide

Vibration Isolation Systems

Optical Tables

Darkrooms/ Dark Boxes

Optical Bases

Posts

Brackets

Adapters

Tools